


Movella Carbon Explosionszeichnung



NO.	Part NO.	Description	Qty
1	HF-B00HR2020001	Right Handle Assembly	1
2	HF-B00HR2020002	Right Frame Assembly	1
3	HF-B00HR2020003	Crossbar Assembly	1
4	HF-B00HR2020004	Storage Bag Assembly	1
5	HF-B00HR2020005	Left Handle Assembly	1
6	HF-ZPSXHR20220A	HR202-20A Brake Cable Assembly (Version B)	2
7	HF-B00HR2020006001	Double-Lock Liner Tube Assembly (Right)	1
8	HF-B00HR2020018001	Double-Lock Liner Tube Assembly (Left)	1
9	HF-B00HR2020007	Backrest Strap Assembly	1
10	HF-B00HR2020008	Side Bracket Assembly	1
11	HF-B00HR2020009	Rear Wheel Assembly	2
12	HF-B00HR2020010	Left Frame Assembly	1
13	HF-B00HR2020011	Front Fork Assembly	2